


Scaling smalldata_tools for LCLS-II

Vincent Esposito and Silke Nelson, LCLS Experiment Control System

Oct 21, 2025

What is smalldata_tools?

From data acquisition to user files



- Configurable pre-processing tool built on top of psana. Two approaches:
 - Event-based reduction
 - Binned reduction
- Output files (hdf5) suited for Jupyter Notebooks analysis / exploration
- Abstract event loop, calibrations, scaling from users
- Used at most instruments at LCLS
- **Users can focus on their science rather than the computing and data processing**

Approach 1: event-based reduction

Keep all events, but make them small


- Area and waveform detectors:
 - Region of interest, projection
 - Azimuthal integration
 - Droplet / photonization
 - Autocorrelation
 - ... (easily extendable)
- Default set of scalar detector

Simple pipeline example: droplet + greedy-guess photon finding algorithms.

```
# Instantiate detector
det = DetObject(detname, run, **det_kwargs)

# Instantiate the different analysis functions
roi_func = ROIFunc(**roi_args)
droplet_func = dropletFunc(**droplet_args)
photon_func = droplet2Photons(**photon_args)
sparsify = sparsifyFunc()

# Build a simple pipeline and add it to the detector
photon_func.addFunc(sparsify)
droplet_func.addFunc(photon_func)
roi_func.addFunc(droplet_func)
det.addFunc(roi_func)
```



Approach 2: cube

Aggregate events, but keep full detectors

- Event aggregation:
 - Step scans
 - User-defined multi-dimensional binning
- Event screening and filtering
 - Laser on / off
 - Discard / filter events

```
def screener(run):
    ddets = hutch_default.rlxDetectors(run)
    lightstatus_det = [d for d in ddets if d.name == "lightStatus"][0]

    # Individual filters
    laser_on = BoolFilter(
        lightstatus_det, "laser", expected_state=True, label="laser_on"
    )

    laser_off = BoolFilter(
        lightstatus_det, "laser", expected_state=False, label="laser_off"
    )

    xray_on = BoolFilter(
        lightstatus_det,
        "xray",
        expected_state=True,
    )

    xray_off = BoolFilter(
        lightstatus_det, "xray", expected_state=False, label="dropped_shots"
    )

    # Combine individual filters into composite filters for different states
    screener_on = CompositeFilter([laser_on, xray_on])
    screener_off = CompositeFilter([laser_off, xray_on])


    # Combine all screeners into an OR screener to be run by the cube
    event_screener = CompositeFilter(
        [screener_on, screener_off, xray_off], require_all=False
    )

    return event_screener
```

Scaling: psana







MPI communicators and groups in psana

- LCLS-1 psana: embarrassingly parallel
- LCLS-2 psana: complex parallelization with multiple stages:
 - SMD0: distribute block of events
 - EB (event builder)
 - BD (big data)
 - SRV: writers
- Knobs to tune each stage of the parallel processing















Run with slurm: a simple interface

The Automated Run Processing (ARP)

| Logbook for rix101265125 | | | | | | |
|--------------------------|------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|----------|--------------|----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Name | Executable | Parameters | Location | Trigger | As user | |
| cube | /sdf/data/lcls/ds/rix/rix101265125/results/smalldata_tools/arp_scripts/lets_cube_lcl2.sh | --config crix | S3DF | MANUAL | jjoshi |   |
| run_summary | /sdf/data/lcls/ds/rix/rix101265125/results/smalldata_tools/arp_scripts/submit_plots.sh | --queue milano | S3DF | MANUAL | dgarratt |   |
| smd | /sdf/data/lcls/ds/rix/rix101265125/results/smalldata_tools/arp_scripts/submit_smd2.sh | --partition milano --postRunTable --nodes 3 --wait --config crix --gather_interval 50 | S3DF | START_OF_RUN | dgarratt |   |


Submit slurm jobs from the elog web interface. No CLI needed.

Processed files available minutes after data run is over

| Logbook for rix101265125 | | | | | | |
|--------------------------|----------------------------------|--------|----------|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------|
| Run | Job | Status | Job ID | Submit Time | Actions | Report |
| 414 | <input type="text" value="smd"/> | DONE | 13364601 | Oct/9/2025 08:39:10 |    | • Last Event: ~ 360 cores * 464 evts |
| 413 | <input type="text" value="smd"/> | DONE | 13207636 | Oct/6/2025 05:54:05 |    | • Last Event: ~ 360 cores * 7386 evts |
| 412 | <input type="text" value="smd"/> | DONE | 13207560 | Oct/6/2025 05:48:56 |    | • Last Event: ~ 360 cores * 11742 evts |
| 411 | <input type="text" value="smd"/> | DONE | 13207468 | Oct/6/2025 05:43:43 |    | • Last Event: ~ 360 cores * 19120 evts |

Challenges and outlook

- LCLS welcomes users with a range of computing skills
- Build trust: more complex systems become black-boxes to the end users
- Balance between stability / continuity and versatility



An aerial photograph of the SLAC National Accelerator Laboratory. A long, straight building with a grey roof and white rectangular panels runs through the center of the image. The building is flanked by dense green trees and a road. In the background, there are rolling hills and a cloudy sky with some sunlight breaking through. The text "Thank you for your attention" is overlaid in white, sans-serif font across the upper middle of the image.

Thank you for your attention